Acetal Discs

Disco de Acetal

Disque de Acetal

Distributed by (in US only)  
HENRY SCHEIN Inc.  
135 DURYEA ROAD  
Melville, NY 11747 USA

HENRY SCHEIN U.K. HOLDINGS LTD.  
MEDICARE HOUSE, CENTURION CLOSE  
GILLINGHAM BUSINESS PARK  
GILLINGHAM, ME8 0SB U.K.  
www.zirlux.com  
MADE IN USA/Fabricado en los EU/US  
Fabriqué aux Etats-Unis  
Rev: 2015/07

## Zirlux Acetal Instructions for Use

### SECTION 1: SYSTEM OVERVIEW

Zirlux Acetal is a pure, highly crystalline acetal copolymer resin possessing high tensile and flexural strength, fatigue resistance, and hardness. These properties, when combined with its very low moisture absorption, make the acetal material ideally suited for milling a variety of removable dental appliances, including partial denture frameworks, provisional crowns and bridges, and bite splints.

### SECTION 2: PRODUCT COMPOSITION

Ingredient	Approx % weight
Polyoxymethylene Copolymer	99 to 100
Pigments	< 1.0

### SECTION 3: TECHNICAL DATA

Property	Method	Unit	Value
Density	ASTM D792	g/cm <sup>3</sup>	1.41
Tensile Strength (New)	ASTM D638	MPa	61.00
Tensile Modulus	ASTM D638	MPa	2820.00
Flexural Strength (New)	ASTM D790	MPa	90.00
Flexural Modulus	ASTM D790	MPa	2620.00
Moisture Absorption	ASTM D570	wt%	0.22
Izod Impact, Notched	ASTM D256	ft-lb/in (J/m)	1.00

### SECTION 4: INDICATIONS

Removable prosthetic dental appliances
Partial denture clasp
Partial denture framework
Provisional crown & bridge
Bite splint

### SECTION 5: CONTRAINDICATIONS

Contraindications are all applications not included in SECTION 4 above.

### SECTION 6: PREPARATION GUIDELINES

General preparation guidelines for all removable appliances should be followed.

### SECTION 7: DESIGN OF RESTORATION

Partial denture framework (minimum thickness)

- Major connector: 0.7 mm
- Saddle area: 0.5 mm
- Clasp on anterior tooth
- 2.0 mm where clasp joins framework, thinning to 0.7 mm at tip

Clasp on posterior tooth

- 2.5 mm where clasp joins framework, thinning to 0.7 mm at tip

Temporary Crown & Bridge

- 0.5 mm minimum thickness
- 0.03 mm cement gap

### SECTION 8: MILLING

1. Mount the disc according to the CAM-System instructions. Set the enlargement-factor into the unit at a ratio of 1:1.
2. Start the milling unit.
3. Select appropriate milling strategy depending on the appliance to be fabricated.
4. Once the milling operation has completed, remove the disc containing the milled restoration(s).

### SECTION 9: POST MILLING

1. Remove the restoration(s) from the disc using a separating disc or suitable acrylic bur.
2. To polish the restoration, load a new polishing wheel or buff with polish suitable for thermoplastic materials and buff all areas of the appliance until all scratches are removed and high shine is achieved.  
Notice: All Zirlux Acetal discs are solely for use by or on the order of a dental professional. They are not for use by the general public or over-the-counter.

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