Instructions for use for Zirlux FC2 Discs

1. PRODUCT DESCRIPTION:

Zirlux FC2 blanks are made of pre-sintered yttria stabilized tetragonal zirconia polycrystalline ceramics (YTZP) for use as CAD/CAM milling blanks. YTZP or 3Y-TZ is a commercially pure zirconia (ZrO_2) comprising about 3 mol % of Y_2O_3 to stabilize the tetragonal phase, the strongest and toughest form of zirconia. Zirlux FC2 blanks are pre-shaded and in a "chalk-like" state and can be milled using any compatible CAD/CAM machine. It is necessary to mill the blanks with an appropriate enlargement factor to account for the shrinkage that occurs during full sintering. The blanks are manufactured to have a well-controlled enlargement factor of approximately 1.22 (22%), which ensures that a good fit is achieved after full sintering. Each blank is labeled with its individual, precisely determined, enlargement factor.

After full sintering, the material densifies into a high strength dental ceramic listed as Type 2 Class 6 according to ISO6872:2008

Indications

- Full Contour anterior and posterior crowns, bridges*, inlays & onlays
- Single tooth and bridge* frameworks in the anterior and posterior

Contraindications

- Improper framework design

COMPOSITION (in wt. %)

ZrO_2	(+ HfO_2):>94%
Y_2O_3	$5.35\% \pm 0.20$
Al_2O_3	<0.1%
HfO_2	Typically <3.0%

TECHNICAL DATA †

Coefficient of Thermal Expansion (25-500°C):	$10.6 \times 10^{-6}/^\circ\text{C}$
Flexural Strength:	1100 MPa
Vickers Hardness (HV10):	$\geq 1,250$
Chemical Solubility:	$\leq 1 \mu\text{g}/\text{cm}^2$

†As per ISO 6872:2008

2. ZIRCONIA FRAMEWORK DESIGN:

When designing a full contour zirconia restoration, it is important that minimum thickness requirements are met (see Table 1). 1.0 mm of occlusal and axial thickness is considered ideal for full contour zirconia restorations.

Table 1: Minimum thickness for full contour zirconia, copings, and frameworks

Anterior	Crowns		Bridges	
	Single-unit	3-unit	≥ 4 units with max 2 connected pontics*	Cantilever with 1 pontic
Axial:	0.5 mm	0.5 mm	0.7 mm	0.7 mm
Incisal/Occlusal:	0.5 mm	0.7 mm	1.0 mm	1.0 mm
Connector Dimensions				
n/a				
Posterior				
Axial:	0.5 mm	0.5 mm	0.7 mm	0.7 mm
Incisal/Occlusal:	0.5 mm	0.7 mm	1.0 mm	1.0 mm
Connector Dimensions:	n/a	9 mm ²	12 mm ²	12 mm ²

*In Canada, zirconia full contour bridges and zirconia bridge frameworks are indicated for 6 units or less with a maximum 2 connected pontics.

3. SHADE SELECTION:

Final Shade	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Disc Shade	U2	U3	U3	U3	U5	U2	U2	U4	U4	U2	U5	U5	U5	U2	U3	U5

Final Shade	SB	A00	B00
Disc Shade	U1	U1	U1

4. MILLING:

Milling of the blanks is done with an enlargement of the zirconia factor of proximately 1.22 (or 22%) relative to the final desired dimensions. This compensates for the shrinkage that occurs during full sintering. Each Zirlux FC2 blank is labeled with its individual and unique Enlargement Factor (EF), as well as the Linear Shrinkage (LS, %) which is calculated from the EF. Please refer to the operating instructions of the corresponding CAD/CAM system for detailed milling instructions and the required inputs. It is strongly recommended to use a 0.5 mm diameter or finer tool for finishing occlusal surfaces to achieve accurate rendering of tooth anatomy and smoother contacts with opposing dentition.

5. FINISHING AND PREPARING FOR SINTERING:

It is recommended that all adjustments to zirconia restorations be completed prior to sintering. Only use suitable grinding instruments, low speed, and light pressure, to avoid flaking, chipping and cracking of the restoration.

Steps:

- Carefully separate the milled restoration from the starting blank using a fine diamond disk and smooth out the attachment area using a suitable grinding instrument.
- Make adjustments to the milled restoration as necessary.
- After finishing, clean the milled restoration using compressed air or water as necessary.
- If the restoration is moist or wet it should be dried prior to full sintering.

6. SINTERING:

The sintering temperature for the zirconia is $1500^{\circ}\text{C} \pm 25^{\circ}\text{C}$. Higher temperature (up to 1525°C) is required for faster sintering cycles. Restorations can be sintered in any dental zirconia sintering furnace provided it can run the recommended sintering cycle below. Follow the sintering instructions of the corresponding furnace with regard to placement of frameworks and furnace accessories such as crucibles, sintering beads, sintering trays, etc. If using a crucible or dish, it is recommended to sinter Zirlux restorations uncovered.

CAUTION: Do not sinter Zirlux restorations along with any restorations stained with coloring liquids. Shade contamination may occur.

Sintering Cycle*

Heating Rate 1:	10°C/m
Temperature 1:	700°C
Hold Time 1:	0:00
Heating Rate 2:	2-5°C/min **
Temperature 2:	1500 °C
Hold Time 2:	2 hours
Cooling:	Furnace/Natural Cool***

*For dry milled restorations. Include a drying step for wet milled restorations.

** Larger mass restorations require a slower heating rate.

***Do not remove restorations from the furnace until below 500°C

If using a sintering furnace with molybdenum disilicide (MoSi_2) heating elements a Cleaning Cycle should be run if: (i) it is the first use of a new heating element(s), or (ii) sintered frameworks have come out discolored (yellowish). The recommended Cleaning Cycle is: heat at $10^{\circ}\text{C}/\text{min}$. to 1550°C , hold for 4 hours, followed by natural cool. The furnace should be run empty for this cycle. The effect of the Cleaning Cycle is to generate a protective silica (SiO_2) coating on the heating elements, as is necessary to prevent discoloration of restorations.

7. POST-SINTERING PREPARATION:

If adjustments to the zirconia are necessary, it is recommended to use only grinding instruments indicated for zirconia. When grinding, use little or no pressure to reduce heat which can cause fractures in the framework. Ensure that the minimum thickness requirements and connector dimensions are maintained. See Tables 1 and 2. It is critical that the connector areas are not touched by a separating disc or grinding tool. This will compromise the strength of the zirconia framework.

Always use water during grinding to keep restorations cool. Before staining, the surface of the zirconia restoration must be sandblasted with alumina not coarser than 50 psi. After sandblasting, steam clean or ultrasonically clean the framework for 15 min. It is very important that the surface of the zirconia is free from dirt, milling dust/residue, and oily/greasy elements. When applying stains to the surface of zirconia, use only those intended for this purpose.

8. STAIN AND GLAZE OR LAYERING:

Refer to the respective Zirlux FC2 Characterization and Zirlux LC instructions for use for post sintering requirements, techniques and firing parameters.

9. POLISHING FULL CONTOUR RESTORATIONS:

It is very important that the occlusal surface of the full contour Zirlux FC2 restoration has a smooth, glazed or polished surface. An abrasive surface can wear the patient's opposing dentition. If any anatomical adjustments are made on zirconia by the laboratory or dentist, the surface must be either re-glazed or polished using polishing cups and pastes developed to finish zirconia to maximum surface smoothness. Before staining, it is not necessary to polish the surfaces which do not come into occlusion with opposing dentition.

**Please note that it is possible to apply Zirlux FC2 Glaze paste over a polished zirconia surface.*

10. STORAGE:

Store in a cool, dry place.

11. WARNING AND SAFETY INSTRUCTIONS:

Do not inhale grinding dust. Do not get grinding dust into eyes. A dust mask and protective eyewear should be worn when processing zirconia restorations.

Manufactured for:

 Ardent, Inc.
175 Pineview Dr.
Amherst, NY 14228
800-922-5520

 AB Ardent
Generatorgatan 8
SE-195 60 Arlandastad
SWEDEN

  R ONLY
For Dental Use/
Made in USA

©2012 Ardent, Inc. All rights reserved.
Zirlux is a registered trademark of Henry Schein, Inc.
REF 640017 Rev. 1 2012-11
For more information, contact Zahn Dental at 1-800-496-9500

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©2012 Ardent, Inc. Alle rechten voorbehouden.
Zirlux is een geregistreerd handelsmerk van Henry Schein, Inc..